



AM-350 Treatment Table

- | | | | |
|--|---|---|--|
| <ul style="list-style-type: none"> • 3-Section Treatment Table • 76"L x 27"W | <ul style="list-style-type: none"> • Min. Height: 18" • Max Height: 37" • 1 1/2" Firm Density Foam | <ul style="list-style-type: none"> • Head: 15" x 27" with contoured face opening • Center: 22" x 27" • Foot: 39" x 27" | <ul style="list-style-type: none"> • 400lbs. Capacity • Power-assisted Caster System |
|--|---|---|--|



AM-353 Treatment Table

- | | | | |
|--|---|---|---|
| <ul style="list-style-type: none"> • 3-Section Treatment Table • 76"L x 27"W | <ul style="list-style-type: none"> • Min. Height: 17" • Max Height: 36" • 1 1/2" Firm Density Foam | <ul style="list-style-type: none"> • Head: 15" x 27" with contoured face opening • Center: 22" x 27" • Foot: 39" x 27" | <ul style="list-style-type: none"> • 400lbs. Capacity • Non-skid Adjustable Glides (No Casters) |
|--|---|---|---|