

Machine Stand Accommodates All Traction Machines



- Friction Free Lumbar Separates 5" And Rolls On 4 Steel Ball Roller Bearings.
- The Rolling Section Can Be Locked In Any Position With A Positive Gear Rack System.

AM-420 Traction Table

- | | | | |
|---|---|--|--|
| <ul style="list-style-type: none"> • 2-Section Traction Table • 76"L x 27"W • 93" Overall Length | <ul style="list-style-type: none"> • Fixed Height: 30" • Shelf: 12" x 66" (upholstered) • 1 1/2" Firm Density Foam | <ul style="list-style-type: none"> • Fixed Section: 34" x 27" • Rolling Section: 42" x 27" | <ul style="list-style-type: none"> • 400lbs. Capacity • Non-skid Adjustable Glides |
|---|---|--|--|



AM-850 Flexion Stool

- | | | |
|-------------------------|-------------------------------------|----------------|
| 1-Section Flexion Stool | Min. Height: 12"
Max Height: 18" | Pad: 12" x 19" |
|-------------------------|-------------------------------------|----------------|

AM-855 Deluxe Flexion Stool

- | | | |
|--|-------------------------------------|------------------------|
| 2-Section Flexion Stool with rubber cups | Min. Height: 13"
Max Height: 19" | Two Leg Pads: 9" x 12" |
|--|-------------------------------------|------------------------|

