



AM-500 Treatment Table

- | | | | |
|--|---|---|--|
| <ul style="list-style-type: none"> • 3-Section Treatment Table • 3-Piece Head • 76"L x 27"W | <ul style="list-style-type: none"> • Min. Height: 18" • Max Height: 37" • 1 1/2" Firm Density Foam | <ul style="list-style-type: none"> • Head: 15" x 27" with contoured face opening • Arms: 12.5" x 5.5" • Center: 22" x 27" • Foot: 39" x 27" | <ul style="list-style-type: none"> • 400lbs. Capacity • Power-assisted Caster System |
|--|---|---|--|



AM-550 Treatment Table

- | | | | |
|--|---|---|--|
| <ul style="list-style-type: none"> • 3-Section Treatment Table • 3-Piece Head • 76"L x 27"W | <ul style="list-style-type: none"> • Min. Height: 18" • Max Height: 37" • 1 1/2" Firm Density Foam | <ul style="list-style-type: none"> • Head: 15" x 27" with contoured face opening • Arms: 12.5" x 5.5" • Center: 22" x 27" • Foot: 39" x 27" | <ul style="list-style-type: none"> • 400lbs. Capacity • Power-assisted Caster System |
|--|---|---|--|