



Adjustable Back
ELEVATES 80°

Matching Wood Shelf
20" x 66"

AM-616 Treatment Table with adjustable back

- 2-Section Treatment Table
- 78"L x 30"W
- Fixed Height: 30"
- 1 3/4" Firm Density Foam
- Backrest: 30" x 31"
- Body: 30" x 47"
- 500lbs. Capacity



Adjustable Back
ELEVATES 80°

Center Leg with Adjustable Glide
for Extra Stability

AM-618 Treatment Table with cabinet and adjustable back

- 2-Section Treatment Table
- 78"L x 30"W
- Fixed Height: 30"
- 1 3/4" Firm Density Foam
- Backrest: 30" x 31"
- Body: 30" x 47"
- 500lbs. Capacity