



AM-SX4000 Treatment Table

- | | | | |
|---|---|--|---|
| <ul style="list-style-type: none"> • 4-Section Treatment Table • 72"L x 27"W • 99" Overall | <ul style="list-style-type: none"> • Min. Height: 18" • Max Height: 37" • 1 1/2" Firm Density Foam | <ul style="list-style-type: none"> • Head: 15" x 27" with contoured face opening • Thoracic: 18" x 27" • Center: 12" x 27" • Foot: 31" x 27" | <ul style="list-style-type: none"> • 400lbs. Capacity • 4 Total Lock Swivel Casters |
|---|---|--|---|



AM-SX1060 Changing Table

- | | | |
|---|---|---|
| <ul style="list-style-type: none"> • 1-Section Changing Table with Side Rails • 60"L x 27"W | <ul style="list-style-type: none"> • Min. Height: 18" • Max Height: 37" • 1 1/2" Firm Density Foam | <ul style="list-style-type: none"> • 400lbs. Capacity • 4 Total Lock Swivel Casters |
|---|---|---|

Also available in Model: AM-SX1072 measuring 72"L x 25"W