NAME:	PHONE:	
SCHOOL/ORGANIZATION:	E-MAIL:	
DEALER/REP:	SALES REP:	QUOTE:



Designed to support the busy practitioner and multiple patients, solid hardwood construction means our Classic Wood Mat Table is built for strength and stability during practitioner-assisted exercise and other treatments. Our 6 Leg Mat table comes standard with solid hardwood on the base frame, six legs and supports with a thick plywood top construction. The mat table also features a layer of comfortable 1.5" high density foam with softly rounded corners and triple-bolted leg assembly for extra durability.

STANDARD FEATURES:

- Load Capacity: 700 lbs.
- Solid hardwood frame and 6 legs
- 3/4" Plywood top
- 1.5" High density foam cushion
- Rounded corners
- Triple-bolted leg assembly
- ADA Compliant

b Leg Classic Wood MAT TABLES ple sic ing Dur on look of leded MAT-4884-6L

STANDARD TABLES:

MAT-4884-6L	MAT-6084-6L	MAT-7296-6L
48"W x 84"L x 19"H	60"W x 84"L x 19"H	72"W x 96"L x 19"H
Oty	Oty	Oty

STANDARD WOOD CHOICES:

	Oak	



WOOD STAIN OPTION:

MAT-Stain

STAIN CHOICES:

Warm Coppe
Deep Garnet
Fossil Gray
Tungsten

Rustic Fall
Warm Caramel
Black Beauty
Custom Wood Stain Color

OPTION:

LOGO-BB - Table Sling



STANDARD UPHOLSTERY CHOICES:

Royal Blue		American Beauty
Imperial Blue		Lark
Blue Ridge		Burgundy
Deep Violet		Dove
Forest		Charcoal
Black		Mocha
Mandarin Orang	ge	