NAME:	PHONE:	
SCHOOL/ORGANIZATION:	E-MAIL:	
DEALER/REP:	SALES REP:	QUOTE:



Beautifully portable, the Basic Pro portable chiropractic table is still tough enough to withstand the requirements of the most demanding use. It features both thoracic and pelvic drops, narrow shoulder design, five angle adjustable headpiece, adjustable ankle extension, 2" multi-layer firm foam cushion, adjustable height and folding legs. The Basic Pro portable chiropractic table comes standard with paper holder, hand rests, and carry case.

STANDARD FEATURES:

- Load Capacity: 400 lbs.
- Table weight: 46 lbs.
- Thoracic & pelvic drops
- Narrow shoulder design
- 5 Angle adjustable headpiece
- Adjustable ankle extension
- Adjustable height, folding legs
- Height range: 18.5-22.4"
- Folded table: 21.5"W x 37"L x 8"H
- 2" Multi-layer foam cushion
- Black upholstery
- Hand rests
- Paper holder
- Black carry case

STANDARD TABLE:

BASIC PRO 21.5"W x 70"L x 18.5-22.4"H

Basic Pro Portable CHIROPRACTIC TABLES









Adjustable Height Folding Legs





Height Range: 18.5-22.4"

Adjustable Headpiece