NAME:	PHONE:	
SCHOOL/ORGANIZATION:	E-MAIL:	
DEALER/REP:	SALES REP:	QUOTE:



Designed to support the busy practitioner and multiple patients, solid hardwood construction means our Classic Wood Mat Table is built for strength and stability during practitioner-assisted exercise and other treatments. Our 4 Leg Mat table comes standard with solid hardwood on the base frame, four legs and supports with a thick plywood top construction. The mat table also features a layer of comfortable 1.5" high density foam with softly rounded corners and triple-bolted leg assembly for extra durability.

STANDARD FEATURES:

- Load Capacity: 450 lbs.
- Solid hardwood frame and 4 legs
- 3/4" Plywood top
- 1.5" High density foam cushion
- Rounded corners
- Triple-bolted leg assembly
- ADA Compliant

Classic Wood MAT TABLES MAT-4872-41

STANDARD TABLES:

MAT-4872-4L	MAT-4884-4L	MAT-6084-4L
48"W x 72"L x 19"H	48"W x 84"L x 19"H	60"W x 84"L x 19"H
Qty	Qty	Qty

STANDARD WOOD CHOICES:

	Oak		Maple	

WOOD STAIN OPTION:

MAT-Stain

STAIN CHOICES:

Warm Copper	Rustic Fall
Deep Garnet	Warm Caramel
Fossil Gray	Black Beauty
Tungsten	Custom Wood Stain Colo

OPTION:

LOGO-BB - Table Sling



STANDARD UPHOLSTERY CHOICES:

Royal Blue	American Beauty
Imperial Blue	Lark
Blue Ridge	Burgundy
Deep Violet	Dove
Forest	Charcoal
Black	Mocha
Mandarin Orange	