NAME:	PHONE:	
SCHOOL/ORGANIZATION:	E-MAIL:	
DEALER/REP:	SALES REP:	QUOTE:



The Convertible Staircase features straight and corner configurations with a 36"W x 36"D x 24"H platform, (4) 36"W x 9"D x 6"H steps and (6) 36"W x 9"D x 4"H steps. The uprights are constructed of 2" steel and the handrails are made from solid hardwood. Anti-slip safety treads are installed on the platform and all steps for increased safety. ADA compliant staircases are also available that include 11"D steps and handrails with 12" of overhang.

STANDARD FEATURES:

- Load Capacity: 350 lbs.
- Stairs can be configured as a straight or corner unit
- Platform is 36"W x 36"D x 24"H
- (4) 36"W x 9"D x 6"H steps and (6) 36"W x 9"D x 4"H steps
- Handrails are made from solid hardwood (ship ready to assemble)
- Uprights are made from 2" steel with a black powder coat finish
- Stairs are available in wood or laminate finishes
- Anti-slip safety tread on platform and all steps for increased safety
- ADA compliant staircases are available with 11"D steps and handrails with 12" of overhang
- 9' ceiling height required



STAIRCASE CONFIGURATIONS:

	PTST-005-WHR
--	--------------

Platform: 36"W x 36"D x 24"H Stairs: 109"L x 36"W x 24"H

9"D Steps

Overall including Handrails (straight): 109"L x 42"W x 59"H Overall including Handrails (corner): 78"L x 60"W x 59"H

Qty _____

☐ PTST-005-WHR-ADA

Platform: 36"W x 36"D x 24"H Stairs: 125"L x 36"W x 24"H

11"D Steps and 12" Overhang on Handrails

Overall including Handrails (straight): 157"L x 42"W x 59"H Overall including Handrails (corner): 102"L x 84"W x 59"H

Qty _____

STANDARD WOOD CHOICES:



STANDARD LAMINATE CHOICES:

