NAME:	PHONE:	
SCHOOL/ORGANIZATION:	E-MAIL:	
DEALER/REP:	SALES REP:	_QUOTE:



The ET Treatment Cabinet sets the standard in functionality, ease of use, storage and convenience. With an easy-to-operate adjustable height ranging from 29" to 44", and the numerous cabinet and top options, this cabinet serves most therapists' needs. The cushion can be upgraded to include 70° lift back and 30° split leg options on the cushion top, and you can upgrade the cabinet configuration to include storage on one or two sides. Enhance your room by adding one of our custom branding options on the cushion or cabinet.

STANDARD FEATURES:

- Load Capacity: 500 lbs.
- 30"W x 83"L x 29-44"H
- Durable metal cabinet with black powder coat
- Leveling feet on cabinet base
- 2" casters for easy movement
- (1) Mounted hand control
- Quiet and smooth lift columns elevate athletes with ease
- 2" High density foam cushion
- Rounded corners

Qty_

- Door locks included on ET-101, ET-102, ET-106, and ET-107 models
- Fixed height (32") treatment cabinets also available



SIANDARD	ELECTRIC	CABINE I:	
----------	----------	-----------	--

	· · · · · · · · · · · · · · · · · · ·	
ET-100 30"W x 83"L x 29-44"H Open Shelf	■ ET-101 30"W x 83"L x 29-44"H Double Doors on One Side ■ Right Side ■ Left Side	☐ ET-105 30"W x 83"L x 32"H Open Shelf
Oty	Oty	Oty
ET-102 30"W x 83"L x 29-44"H Double Doors on Both Sides	■ ET-103 30"W x 83"L x 29-44"H Sliding Doors on One Side □ Right Side □ Left Side	ET-107 30"W x 83"L x 32"H Double Doors on Both Sides
Oty	Oty	Oty
ET-104 30"W x 83"L x 29-44"H Sliding Doors on Both Sides		ET-109 30"W x 83"L x 32"H Sliding Doors on Both Sides

TANDARD FIXED HEIGHT CABINET:			
ET-105 30"W x 83"L x 32"H Open Shelf	■ ET-106 30"W x 83"L x 32"H Double Doors on One Side ■ Right Side ■ Left Side		
Qty	Oty		
ET-107 30"W x 83"L x 32"H Double Doors on Both Sides	■ ET-108 30"W x 83"L x 32"H Sliding Doors on One Side ■ Right Side ■ Left Side		
Qty	Oty		
ET-109 30"W x 83"L x 32"H			

Qty _____

5



ET TREATMENT CABINET

OPTIONS:

- ET-201 70° Lift Back Cushion
- ET-202 Flat, 30° Split Leg Cushion
- **ET-203** 70° Lift Back Cushion, 30° Split Leg Cushion

STANDARD UPHOLSTERY CHOICES:

Royal Blue	American Beauty
Imperial Blue	Lark
Blue Ridge	Burgundy
Deep Violet	Dove
Forest	Charcoal
Black	Mocha
Mandarin Orange	

STANDARD FRAME COLORS:

Black

CUSTOM FRAME COLORS:

White	Green
Light Gray	Forest Green
Charcoal	Purple
Royal Blue	Yellow
Navy Blue	Orange
Red	Brown
Burgundy	Bright Blue

ET-103 with Lift Back, Split Leg, Color Edge Print Logo, Custom Frame Color and Large Logo Decal (Doors on Left Side)



CUSTOM BRANDING:

Decal Custom Small Small Logo Decal (6"W x 6"H)



Decal Custom Large Large Logo Decal (12"W x 12"H)



CUSTOM CUSHION BRANDING:

Logo Slip Cover (Lift Back Only) (30"W x 12"H)



Deboss Logo (12"W x 12"H) DEBPLT12 Deboss Plate



LOGO CP Color Edge Print (12"W x 12"H)



COLOR EDGE PRINT UPHOLSTERY CHOICES:



**Color Edge Print upholstery colors do not match Standard Upholstery colors.