NAME:	PHONE:	
SCHOOL/ORGANIZATION:	E-MAIL:	
DEALER/REP:	SALES REP:	QUOTE:



Our space-saving Portable Wood Mat Table provides an ideal surface for stretching, therapeutic exercise or massage, yet it can be quickly and easily stored out of the way for more room when you have limited space. This rugged, 800-pound capacity table features premium materials and unmatched craftsmanship. Built using only the finest components such as eastern hard rock maple, American poplar, extra rigid birch and aircraft-quality cables, it is able to withstand the most vigorous work. The Cable-Lock™ system provides extraordinary strength and easy setup, and the unique Inverted Truss design gives extra support in the middle of the table, ensuring a solution that is both foldable and ruggedly stable. The 1" extra firm foam configuration is ideal for a variety of treatment protocols, and the adjustable height legs feature two stabilizing knobs to eliminate table movement.

STANDARD FEATURES:

- Load Capacity: 800 lbs.
- Table weight: 41 lbs.
- Low height range: 17-25"
- Solid extra rigid birch frame with inverted truss
- Cable Lock system
- Full length piano hinge
- Vinyl coated aircraft cables
- Dual knobs
- Reiki end panels
- 9" on center dual-end accessory ports
- 1" Extra firm foam

STANDARD TABLES:

Pmat-3630
36"W x 78"L x 17-25"H
Qty

STANDARD UPHOLSTERY CHOICES:

	Royal Blue		American Beaut
	Imperial Blue		Lark
	Blue Ridge		Burgundy
	Deep Violet		Dove
	Forest		Charcoal
	Mandarin Orange		Mocha
	Black		

Portable Wood MAT TABLES Pmat-3630

OPTIONS:

□ RE - 1.5" Rebond Foam System Medium Density Foam
☐ 19393 - Table Cart
☐ CAXWXL Extra Wide Carry Case
Oty