NAME:	PHONE:	
SCHOOL/ORGANIZATION:	E-MAIL:	
DEALER/REP:	SALES REP:	_QUOTE:



The PHS Medical Portable Treatment Mirror is specifically designed with distortion free glass and is available with a single or double-sided 72" mirror. It can be moved around your facility hassle-free with its 2" swivel casters to accommodate treatment needs.

### **STANDARD FEATURES:**

- Portable full length mirror
- Distortion free glass
- 2" Casters
- Laminate or wood options available

# Portable **MIRRORS** MI-1000 (Back)

# **STANDARD DIMENSIONS:**

MI-100

Oty \_\_\_\_\_

# 24"W x 20"D x 72"H

### **STANDARD WOOD CHOICES:**







## **WOOD STAIN OPTION:**

CAB-Stain

# **STAIN CHOICES:**

Warm Copper Deep Garnet

	Fossil Gray
	Tungsten

Town or I	Rustic Fall
	Warm Caramel

Black Beauty Custom Wood Stain Color

# **OPTIONS:**

MI - Mirror on Back Side

### STANDARD LAMINATE CHOICES:

Oak

Grey Storm Grey

Maple

Harbor Grain

Cherry Merlot Espresso

Black