NAME:	PHONE:		
SCHOOL/ORGANIZATION:	E-MAIL:		
DEALER/REP:	SALES REP:	QUOTE:	



Ease of use, multiple options and a sleek but durable design make the THERA-P line of electric treatment tables the choice of athletic trainers everywhere. Our bariatric model is built strong for larger patients, thanks to precision-cut flat and tube steel. Four top options allow athletic trainers to choose the table that best meets their needs and that ensures patients' most comfortable positioning. It also has a whisper-quiet lift system (with elevating foot control pedal) and easy-to-operate caster system for better ergonomics. Softly rounded corners and a layer of thick 2" foam add the finishing touch.

STANDARD FEATURES:

- Load Capacity: 650 lbs.
- 40"W x 77"L x 17-35"H
- Durable and stylish steel base
- White powder coat base
- Whisper quiet lift system
- Height adjustment: 17-35"
- Single foot control pedal
- 40" Wide cushions in 4 different configuration options
- 2" High density foam cushions
- Rounded corners
- Low profile locking casters
- Head section angle range: +40° to -90°
- Lower extremity angle range: 0 to 65°
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- ADA Compliant

NOTE: International (240V) option available

Thera-PMBariatric ELECTRIC TREATMENT TABLE




TABLE CONFIGURATIONS:



PT1111 1 Section

Qty _____



☐ PT1112 2 Section

Qty _____

PT1113 3 Section

Oty _____

	•	

PT1115 5 Section Qty _____ **OPTIONS:**

PT9025 - 650 lb Lifting Capacity (240V)


■ **\$1136** - Manual Elevating Midsection (Only available on PT1113 & PT1115)


95284 - Dual Foot Control Pedal


\$1100 - Dual Foot Bar Control




TABLE FEATURES:










STANDARD UPHOLSTERY CHOICES:



STANDARD BASE COLORS:

White