

NAME: _____

PHONE: _____

SCHOOL/ORGANIZATION: _____

E-MAIL: _____

DEALER/REP: _____

SALES REP: _____ QUOTE: _____




Ease of use, multiple options and a sleek but durable design make the Thera-P Bobath line of electric treatment tables the choice of therapists everywhere. Four top options allow therapists to choose the table that best meets their needs and that ensures patients' most comfortable positioning. The Thera-P Bobath Electric Treatment Table is constructed of precision-cut flat and tube steel for durability and style. It also has a whisper-quiet lift system with elevating foot pedal.

STANDARD FEATURES:

- Load Capacity: 650 lbs.
- 40"W x 77"L x 17-35"H
- Durable and stylish steel base (PT9024/PT9025)
- White powder coat base
- Whisper quiet lift system with single foot control pedal
- Height adjustment: 17-35"
- 40" Wide cushions in 4 different configuration options
- 2" High density foam cushions with rounded corners
- Low profile locking casters
- Head section angle range: +40° to -90°
- Lower extremity angle range: 0 to 65°
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- ADA Compliant


NOTE: International (240V) option available

Thera-PTM Bobath ELECTRIC TREATMENT TABLE



PT1112

TABLE CONFIGURATIONS:



☐ **PT1111**
1 Section

Qty _____

☐ **PT1112**
2 Section

Qty _____


☐ **PT1113**
3 Section

Qty _____

☐ **PT1115**
5 Section

Qty _____

TABLE FEATURES:




OPTIONS:

- ☐ **PT9025** - 650 lb Lifting Capacity
- ☐ **S1077** - Dual Foot Control Pedal
- ☐ **S1136** - Manual Elevating Midsection (Only available on PT1113 & PT1115)
- ☐ **S1078** - Dual Foot Bar Control



Dual Foot Control Pedal (x2) (S1077)



Dual Foot Bar Control (S1078)



Manual Elevating Midsection (S1135)

STANDARD UPHOLSTERY CHOICES:

- | | |
|--|--|
| ☐ Royal Blue | ☐ American Beauty |
| ☐ Imperial Blue | ☐ Lark |
| ☐ Blue Ridge | ☐ Burgundy |
| ☐ Deep Violet | ☐ Dove |
| ☐ Forest | ☐ Charcoal |
| ☐ Black | ☐ Mocha |
| ☐ Mandarin Orange | ☐ Putty |

STANDARD BASE COLORS:

- ☐ White

CUSTOM BASE COLORS:

- | | |
|-------------------------------------|---------------------------------------|
| ☐ Black | ☐ Green |
| ☐ Light Gray | ☐ Forest Green |
| ☐ Charcoal | ☐ Purple |
| ☐ Royal Blue | ☐ Yellow |
| ☐ Navy Blue | ☐ Orange |
| ☐ Red | ☐ Brown |
| ☐ Burgundy | ☐ Bright Blue |

CUSTOM BASE OPTIONS:

- ☐ **S1178B**
Table Base Graphics
- ☐ **PT9024-CC**
Pick Custom Base Color
- ☐ **S2000B**
Custom Table Base Powder Coat Color and Table Base Graphics

