

Name:	
School/Organization:	
Dealer/Rep:	

Phone:	
E-mail:	
Sales Rep:	Quote:

APOLLO

CLASS IV COLD LASERS

DESKTOP

APOLLO LASER SYSTEM:

- 9"W x 11.25"D x 4.5"H
- 2 Channel desktop unit with choice of probe
- Heavy-duty construction for office and clinic environments
- Probe holders on both sides of the unit that can accommodate two probes and both right and left handed practitioners
- Elevated face for easy viewing and control
- Two timers allow you to treat 2 patients at the same time
- Accessories Included: 12V DC Adapter, Operation Manual, 2 pair of Laser Safety Goggles, 1 Probe Cable, Laser Tutorial DVD, Light & Laser Therapy: Clinical Procedures Book, and 50 Patient Brochures



AP2-DT-5000-C 2-channel desktop unit shown with 5000mW probe

LASER CONFIGURATIONS:

AP2-DT2-5000-C <ul style="list-style-type: none"> • 2 Channel Desktop Unit • 5000mW Cluster Probe 	AP2-DT2-4000-C <ul style="list-style-type: none"> • 2 Channel Desktop Unit • 4000mW Cluster Probe 	AP2-DT2-3000-C <ul style="list-style-type: none"> • 2 Channel Desktop Unit • 3000mW Cluster Probe 	AP2-DT2-500-S <ul style="list-style-type: none"> • 2 Channel Desktop Unit • 500mW Point Probe
QTY:	QTY:	QTY:	QTY:

PORTABLE

APOLLO LASER SYSTEM:

- 5"W x 9"D x 2.5"H
- Portable control unit with choice of probe
- Heavy-duty construction for office and clinic environments
- Weighing only 2.3 pounds, the Apollo laser is the lightest, most powerful portable laser you can buy
- Custom carrying case allows for easy travel to and from office and out in the field
- User replaceable battery
- Accessories Included: Travel bag, 18V DC Adapter, Operation Manual, 2 pair of Laser Safety Goggles, 1 Probe Cable, Laser Tutorial DVD, Light & Laser Therapy: Clinical Procedures Book, and 50 Patient Brochures



AP2-PT-3000-C Portable unit shown with 3000mW probe

LASER CONFIGURATIONS:

AP2-PT-3000-C <ul style="list-style-type: none"> • Portable Control Unit • 3000mW Cluster Probe 	AP2-PT-500-S <ul style="list-style-type: none"> • Portable Control Unit • 500mW Point Probe
QTY:	QTY:

CONTROL UNITS ONLY:

AP2-PT - Apollo AP2-PT Portable Control Unit (Add a portable control unit to your existing desktop unit) Power supply not included.	AP2-DT2 - 2 Channel Desktop Control Unit (Add a desktop control unit to your existing portable unit) Power supply not included.
QTY:	QTY:

PROBES ONLY:

AP2-5000-C - 5000mW Cluster Probe * If probe is to be used with a previous Desktop model there is a software upgrade fee. Compatible with: Desktop Unit	AP2-4000-C - 4000mW Cluster Probe Compatible with: Desktop Unit	AP2-3000-C - 3000mW Cluster Probe Compatible with: Desktop & Portable Units
QTY:	QTY:	QTY:

AP2-500-S - 500mW Point Probe Compatible with: Desktop & Portable Units
QTY:

