Name:	Phone:		
School/Organization:	E-mail:		
Dealer/Rep:	Sales Rep:	Quote:	

## **BASIC PRO** PORTABLE CHIROPRACTIC TABLE

## **STANDARD FEATURES:**

• Load capcity: 400 lbs. • Table weight: 46 lbs. • Height range: 18.5-22.4"

• Folded table: 21.5"W x 37"L x 8"H

• Thoracic & pelvic drops

• Narrow shoulder design

• 5 Angle adjustable headpiece

• Adjustable ankle extension

• Adjustable height, folding legs

• 2" Multi-layer foam cushion

• Black upholstery

• Hand rests

• Paper holder

• Black carry case



## **STANDARD TABLE:**

**BASIC PRO** 21.5"W x 70"L x 18.5-22.4"H

QTY:








Adjustable Height Folding Legs







Adjustable Headpiece

phsmedical solutions.com