

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

# EPD PACKAGE 2

## EPIDURAL POSITIONING DEVICE & STAND

### STANDARD FEATURES:

- Load capacity: 600 lbs.
- 20.5"W x 30"D x 50-60"H (Assembled)
- White powder coat base
- Black upholstery
- Includes 1 Box (50) Sani-Cover® Fitted Disposable Face Rest Covers
- Epidural Positioning Device (EPD)
- 180° Adjustable face rest
- 6 Position adjustable arm rest
- Arm rest size: 16.5"W x 10"D
- Arm rest adjustment: 32-43" from the floor
- Interchangeable torso cushions (1 flat and 1 elevated)
- Torso cushion size: 8"W x 8"H
- Torso cushion adjustment: 5.5-10"
- EPD Stand (LDS)
- Positive locking clamp for EPD
- 7 Position adjustable foot rest
- Foot rest size: 20"W x 12"D
- Foot rest adjustment: 10.5-19.5" from the floor
- Medical grade profile locking casters

EPD Package 2  
(Epidural Positioning  
Device + Stand)  
with flat torso cushion



Elevated torso  
cushion



FRC 30-50



Patents:  
US 6,758,447,B2  
and D455,495

### STANDARD PACKAGE:


<b>EPD Package 2</b> Epidural Positioning Device & EPD Stand
QTY:



### OPTIONS:

<b>FRC 30-50</b> Sani-Cover® Fitted Disposable Face Rest Covers (50 Covers per Box)
QTY:

<b>FRC 30-50-CASE</b> Sani-Cover® Fitted Disposable Face Rest Covers (50 Covers per Box) Case of 18 Boxes (900)
QTY:



EPD - Epidural  
Positioning Device



LDS - EPD Stand



Sani-Covers™ can  
be used on the face  
rest and the arm  
rest of the EPD