



**BASIC**



A Pivotal Health Solutions Brand

[stravasolutions.com](http://stravasolutions.com)

# INTRODUCING A NOT-SO-BASIC SPLINT BATH

**NEW**

*The Strava Basic Splint Bath offers innovative features and modern looks with an economy price*

## Modern Looks

Made of a thermoset material, known for its insulation, flammability, resistance and durability. Features a stainless steel hinged bi-fold lid for easy access from three sides.

## Auto On/Off

Automatically preheats the water to reach a set temperature at a set day and time. Ready when your clinic opens and energy efficient at night when not in use.

## Heats in 25 minutes

Internal heat exchanger is immersed directly in the water basin, ensuring fast heating and uniform heat distribution.

\*\*\*PATENT PENDING\*\*\*







- ✓ Auto On/Off
- ✓ Fast & Uniform Heating
- ✓ User Sets Time & Temperature
- ✓ Audible/Visual Alerts
- ✓ Low Water Sensor
- ✓ Rapid Water Evacuation
- ✓ Bi-Fold Stainless Steel Lid

## DIGITAL TOUCH INTERFACE

**Traditional splint bath controls are antiquated and imprecise**

A full menu of features allow clinicians control Strava's features.



Fully retractable bi-fold lid fits under most cabinets for efficient use of counter top space

