

Name: _____

Phone: _____

School/Organization: _____

E-mail: _____

Dealer/Rep: _____

Sales Rep: _____ Quote: _____

STRAIGHT & PASS THROUGH STAIRCASES

STANDARD FEATURES:

- Load capacity: 350 lbs. (Maple Laminate)
& 450 lbs. (Maple & Birch Wood)
- Anti-microbial and anti-bacterial black powder coat on handrail and uprights
- Uprights are made from 2"W steel
- Stairs are available in wood or laminate finishes
- Anti-slip safety tread on platform and all steps for increased safety
- 9' Ceiling height required
- Staircases ship fully assembled



PTST-001-ALHR



PTST-002-ALHR

STAIRCASE CONFIGURATIONS:

- ☐ PTST-002-ALHR - Straight
Platform: 30"W x 24"D x 24"H
Stairs: 52"L x 30"W x 24"H
(4) 6"H Steps
9"D Treads
Overall Including Handrails: 52"L x 34"W x 61"H

QTY: _____

- ☐ PTST-001-ALHR - Pass Through
Platform: 30"W x 24"D x 12"H
Stairs: 52"L x 30"W x 12"H
(2) 6"H Steps and (3) 4"H Steps
9"D Treads
Overall Including Handrails: 52"L x 34"W x 46"H

QTY: _____

- ☐ PTST-002-ALHR-ADA - Straight ADA
Platform: 36"W x 24"D x 24"H
Stairs: 58"L x 36"W x 24"H
(4) 6"H Steps
11"D Treads
Overall Including Handrails: 58"L x 40"W x 61"H

QTY: _____

- ☐ PTST-001-ALHR-ADA - Pass Through ADA
Platform: 36"W x 24"D x 12"H
Stairs: 58"L x 36"W x 12"H
(2) 6"H Steps and (3) 4"H Steps
11"D Treads
Overall Including Handrails: 58"L x 40"W x 46"H

QTY: _____

STANDARD WOOD CHOICES:


- ☐  Birch
- ☐  Maple

STANDARD LAMINATE CHOICES:

- ☐  Maple



Toll Free: 800.743.7738 | Fax: 605.882.8398
phsmedicalsolutions.com



BUS STEP CONFIGURATIONS:

- ☐ PTST-BUS-STEP
 28.37"W x 9"D x 6"H
 Installed on step is 12"H
 Compatible with the Small Pass Through Staircase and
 Straight Staircase, not the ADA versions
 QTY: _____

STANDARD WOOD CHOICES:

- ☐  Birch
☐  Maple

STANDARD LAMINATE CHOICES:

- ☐  Maple