

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

AM SERIES

27" HI-LO TREATMENT TABLES, NO CASTERS

STANDARD FEATURES:

- Load capacity: 500 lbs
- 27"W x 76"L x 18-37"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1.5" High density foam
- Rounded corners
- Single foot control pedal
- White powder coat base





AM-100

TABLE CONFIGURATIONS:

☐ AM-100 - 1 Section	☐ AM-1227 - 2 Section	☐ AM-353 - 3 Section
QTY	QTY	QTY



TABLE FEATURES:

 Flat Cushion	 Lower Extremity Angle Ranges from 0 to 75°	 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90°
AM-100	AM-1227	AM-353

OPTIONS:

☐ AM-804 - Paper Holder with Cutter (Accommodates 18-21" Paper Rolls)
☐ AM-91 - Add 1" Extra Soft Foam
☐ SPCO - Special Vinyl Color (Provided by Customer)
☐ SPBASE - Special Base Color (Select Custom Base Color)

STANDARD UPHOLSTERY CHOICES:

☐ Imperial Blue	☐ Patina
☐ Blue Ridge	☐ Taupe
☐ Forest Green	☐ Burgundy
☐ Black	☐ Dove Gray
☐ Cappuccino	☐ Greystone

STANDARD BASE COLOR:

☐	Office White
-----------------------	--------------

CUSTOM BASE COLORS:

☐	Black
☐	Steel Gray
☐	Metallic Silver