

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

AM SERIES

27" HI-LO TREATMENT TABLES

STANDARD FEATURES:

- Load capacity: 500 lbs
- 27"W x 76"L x 18-37"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1.5" High density foam
- Rounded corners
- Single foot control pedal
- White powder coat base
- Dual front in-line casters with 3rd wheel drop down for easy maneuverability



AM-350

TABLE CONFIGURATIONS:










☐ AM-150 - 1 Section	☐ AM-200 - 2 Section	☐ AM-227 - 2 Section	☐ AM-250 - 2 Section with 3pc Head Section	☐ AM-350 - 3 Section
QTY	QTY	QTY	QTY	QTY



☐ AM-300 - 3 Section with Elevating Midsection	☐ AM-550 - 5 Section with Tilt Down Armrests	☐ AM-555 - 5 Section with Adjustable Armrests	☐ AM-500 - 5 Section with Elevating Midsection & Tilt Down Armrests	☐ AM-505 - 5 Section with Elevating Midsection & Adjustable Armrests
QTY	QTY	QTY	QTY	QTY



TABLE FEATURES:

 Flat Cushion	 Head Section Angle Ranges from +45° to -90°	 Lower Extremity Angle Ranges from 0 to 78°	 Head Section Angle Ranges from +45° to -90° Armrests Angle Ranges 0° to -45°	 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90°
AM-150	AM-200	AM-227	AM-250	AM-350
 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90° Manual Elevating Midsection	 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90° Armrests Angle Ranges 0° to -45°	 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90° Armrests Vertical Adjustment 8" down	 Armrests Angle Ranges 0° to -45° Head Section Angle Ranges from +45° to -90° Manual Elevating Midsection	 Armrests Vertical Adjustment 8" down Head Section Angle Ranges from +45° to -90° Manual Elevating Midsection
AM-300	AM-550	AM-555	AM-500	AM-505

OPTIONS:

☐ AM-1024 - Add Base with 4 Total Lock Swivel Casters
☐ AM-1060 - Add Base with 4 Casters and 6" Clearance (Adds 4")
☐ AM-804 - Paper Holder with Cutter (Accommodates 18-21" Paper Rolls)
☐ AM-18 - Face Cutout
☐ AM-17 - Face Plug (Must order AM-18 also if ordering AM-150)
☐ AM-91 - Add 1" Extra Soft Foam
☐ AM-807 - Rail Side Fold Down (Pair for 27"W Tables)
☐ AM-808 - Rail Side Fold Down (Pair for 30"W Tables)
☐ AM-875 - Shoulder Mobilization Accessory (27"W, 3 & 5 Section Models Only)
☐ AM-1025 - Adjustable Armrest Package (2 Section Only)
☐ SP30 - 30"W
☐ 19855 - Add-on Belt & Brackets (27" W)
☐ 19859 - Add-on Belt & Brackets (30" W)
☐ SPCO - Special Vinyl Color (Provided by Customer)
☐ SPBASE - Special Base Color (Select Custom Base Color)
☐ 12088 - Side Panel 9 x 33.25 Perforated
☐ 12075 - Side Panel 9 x 33.25 Flat

STANDARD UPHOLSTERY CHOICES:

☐ Imperial Blue	☐ Patina
☐ Blue Ridge	☐ Taupe
☐ Forest Green	☐ Burgundy
☐ Black	☐ Dove Gray
☐ Cappuccino	☐ Greystone

STANDARD BASE COLOR:

☐ Office White

CUSTOM BASE COLORS:

☐ Black
☐ Steel Gray
☐ Metallic Silver