

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

AM-BA SERIES

27" HI-LO TREATMENT TABLES

STANDARD FEATURES:

- Load capacity: 500 lbs
- 27"W x 76"L x 18-37"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- Steel base with dual foot bar control
- 1.5" High density foam
- Rounded corners
- White powder coat base
- (4) Total locking casters
- ADA Compliant



AM-BA350

TABLE CONFIGURATIONS:

☐ AM-BA150 - 1 Section	☐ AM-BA200 - 2 Section	☐ AM-BA227 - 2 Section	☐ AM-BA250 - 2 Section with 3pc Head Section
QTY	QTY	QTY	QTY





☐ AM-BA350 - 3 Section	☐ AM-BA300 - 3 Section with Elevating Midsection	☐ AM-BA550 - 5 Section with Tilt Down Armrests
QTY	QTY	QTY



☐ AM-BA555 - 5 Section with Adjustable Armrests	☐ AM-BA500 - 5 Section with Elevating Midsection & Tilt Down Armrests	☐ AM-BA505 - 5 Section with Elevating Midsection & Adjustable Armrests
QTY	QTY	QTY



TABLE FEATURES:

 Flat Cushion	 Head Section Angle Ranges from +40° to -90°	 Lower Extremity Angle Ranges from 0 to 78°	 Armrests Angle Ranges 0° to -45°	 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90°
AM-BA150	AM-BA200	AM-BA227	AM-BA250	AM-BA350
 Lower Extremity Angle Ranges from 0 to 78° Manual Elevating Midsection Head Section Angle Ranges from +45° to -90°	 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90° Armrests Angle Ranges 0° to -45°	 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90° Armrests Vertical Adjustment 8" down	 Armrests Angle Ranges 0° to -45° Manual Elevating Midsection Head Section Angle Ranges from +45° to -90°	 Armrests Vertical Adjustment 8" down Manual Elevating Midsection Head Section Angle Ranges from +45° to -90°
AM-BA300	AM-BA550	AM-BA555	AM-BA500	AM-BA505

OPTIONS:

☐ AM-804 - Paper Holder with Cutter (Accommodates 18-21" Paper Rolls)
☐ AM-18 - Face Cutout
☐ AM-17 - Face Plug (Must order AM-18 also if ordering AM-BA150)
☐ AM-91 - Add 1" Extra Soft Foam
☐ AM-875 - Shoulder Mobilization Accessory (3 & 5 Section Only)
☐ SP30 - 30"W
☐ 19855 - Add-on Belt & Brackets (27" W)
☐ 19859 - Add-on Belt & Brackets (30" W)
☐ SPCO - Special Vinyl Color (Provided by Customer)
☐ SPBASE - Special Base Color (Select Custom Base Color)
☐ 12088 - Side Panel 9 x 33.25 Perforated
☐ 12075 - Side Panel 9 x 33.25 Flat
☐ 14815 - Nose Hole Diamond

STANDARD UPHOLSTERY CHOICES:

☐ Imperial Blue	☐ Patina
☐ Blue Ridge	☐ Taupe
☐ Forest Green	☐ Burgundy
☐ Black	☐ Dove Gray
☐ Cappuccino	☐ Greystone

STANDARD BASE COLOR:

☐ Office White

CUSTOM BASE COLORS:

☐ Black
☐ Steel Gray
☐ Metallic Silver