

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

AM-BA SERIES

40" BOBATH HI-LO TREATMENT TABLES

STANDARD FEATURES:

- Load capacity: 500 lbs
- 40"W x 76"L x 18-37"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1.5" High density foam
- Rounded corners
- Steel base with dual foot bar control
- White powder coat base
- (4) Total locking casters
- ADA Compliant



AM-BA340

TABLE CONFIGURATIONS:

☐ AM-BA140 - 1 Section	☐ AM-BA240 - 2 Section	☐ AM-BA340 - 3 Section
QTY	QTY	QTY



TABLE FEATURES:

 <p>Flat Cushion</p> <p>AM-BA140</p>	 <p>Head Section Angle Ranges from 0° to 78°</p> <p>AM-BA240</p>	 <p>Head Section Angle Ranges from +45° to -90°</p> <p>Lower Extremity Angle Ranges from 0° to 78°</p> <p>AM-BA340</p>
--	---	---



Armedica®
BY PIVOTAL HEALTH SOLUTIONS

OPTIONS:

☐ AM-17 - Face Plug
☐ AM-91 - Add 1" Extra Soft Foam
☐ 19859 - Add-on Belt & Brackets
☐ SL 72 - Special 72" Length
☐ SPCO - Special Vinyl Color (Provided by Customer)
☐ SPBASE - Special Base Color (Select Custom Base Color)
☐ 12088 - Side Panel 9 x 33.25 Perforated
☐ 12075 - Side Panel 9 x 33.25 Flat

STANDARD UPHOLSTERY CHOICES:

☐  Imperial Blue	☐  Patina
☐  Blue Ridge	☐  Taupe
☐  Forest Green	☐  Burgundy
☐  Black	☐  Dove Gray
☐  Cappuccino	☐  Greystone

STANDARD BASE COLOR:

☐  Office White
--

CUSTOM BASE COLORS:

☐  Black
☐  Steel Gray
☐  Metallic Silver