

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

AM-BAX SERIES

HI-LO TREATMENT TABLES

STANDARD FEATURES:

- Load capacity: 500 lbs
- 27"W x 72 or 76"L x 19-38"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1" High density foam
- Rounded corners
- (4) Total locking casters
- Steel base with dual foot bar control
- White powder coat base



AM-BAX1000

TABLE CONFIGURATIONS:





☐ AM-BAX1000 - 1 Section 27"W x 76"L x 19-38"H	☐ AM-BAX2000 - 2 Section 27"W x 72"L x 19-38"H	☐ AM-BAX2500 - 2 Section with 3pc Head Section 27"W x 72"L x 19-38"H	☐ AM-BAX3500 - 3 Section 27"W x 72"L x 19-38"H
QTY	QTY	QTY	QTY
☐ AM-BAX5400 - 3 Section with Crescent Armrest 27"W x 72"L x 19-38"H	☐ AM-BAX3000 - 3 Section with Electric Elevating Midsection 27"W x 72"L x 19-38"H	☐ AM-BAX5500 - 5 Section with Adjustable Armrests 27"W x 72"L x 19-38"H	☐ AM-BAX5000 - 5 Section with Electric Elevating Midsection & Adjustable Armrests 27"W x 72"L x 19-38"H
QTY	QTY	QTY	QTY



Armedica
BY PIVOTAL HEALTH SOLUTIONS

800.701.5122
armedicamfg.com

TABLE FEATURES:

 Flat Cushion	 Head Section Angle Ranges from +25° to -80°	 Head Section Angle Ranges from +25° to -80°	 Head Section Angle Ranges from +25° to -80°
AM-BAX1000	AM-BAX2000	AM-BAX2500	AM-BAX3500
 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +25° to -80°	 Electric Elevating Midsection Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +25° to -80°	 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +25° to -80° Armrests Vertical Adjustment 8" down	 Armrests Vertical Adjustment 8" down Electric Elevating Midsection Head Section Angle Ranges from +25° to -80°
AM-BAX5400	AM-BAX3000	AM-BAX5500	AM-BAX5000

OPTIONS:

☐ AM-27 - Face Plug
☐ AM-91 - Add 1" Extra Soft Foam
☐ 19855 - Add-on Belt & Brackets
☐ SPCO - Special Vinyl Color (Provided by Customer)
☐ SPBASE - Special Base Color (Select Custom Base Color)

STANDARD BASE COLOR:

☐	Office White
-----------------------	---

CUSTOM BASE COLORS:

☐	Black
☐	Steel Gray
☐	Metallic Silver

STANDARD UPHOLSTERY CHOICES:

☐	Imperial Blue	☐	Patina
☐	Blue Ridge	☐	Taupe
☐	Forest Green	☐	Burgundy
☐	Black	☐	Dove Gray
☐	Cappuccino	☐	Greystone