

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

AM-SP SERIES

HI-LO TREATMENT TABLES

STANDARD FEATURES:

- Load capacity: 400 lbs
- 27"W x 76"L x 24-35.5"H
- One-step caster system
- Single pedestal base
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1.5" High density foam
- Rounded corners
- Single foot control pedal
- White powder coat base






AM-SP100

NOTE: International (240V) option available

TABLE CONFIGURATIONS:

☐ AM-SP100 - 1 Section	☐ AM-SP200 - 2 Section	☐ AM-SP202 - 2 Section with Prenatal Cutout
QTY	QTY	QTY
☐ AM-SP350 - 3 Section	☐ AM-SP300 - 3 Section with Electric Elevating Midsection	☐ AM-SP500 - 5 Section with Electric Elevating Midsection
QTY	QTY	QTY

TABLE FEATURES:

 Flat Cushion	 Head Section Angle Ranges from +45° to -90°	 Head Section Angle Ranges from +45° to -90°
AM-SP100	AM-SP200	AM-SP202
 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90°	 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90° Electric Elevating Midsection	 Lower Extremity Angle Ranges from 0 to 78° Head Section Angle Ranges from +45° to -90° Armrests Range from +45° to -90°
AM-SP350	AM-SP300	AM-SP500

OPTIONS:

☐ AM-SP1024 - Add Base with 4 Total Lock Swivel Casters
☐ AM-18 - Face Cutout
☐ AM-17 - Face Plug (Must order AM-18 also if ordering AM-SP100)
☐ AM-91 - Add 1" Extra Soft Foam
☐ SP30 - 30"W
☐ 19859 - Add-on Belt & Brackets
☐ SPCO - Special Vinyl Color (Provided by Customer)
☐ SPBASE - Special Base Color (Select Custom Base Color)

STANDARD BASE COLOR:

☐	Office White
--------------------------	--------------

CUSTOM BASE COLORS:

☐	Black
☐	Steel Gray
☐	Metallic Silver

STANDARD UPHOLSTERY CHOICES:

☐	Imperial Blue	☐	Patina
☐	Blue Ridge	☐	Taupe
☐	Forest Green	☐	Burgundy
☐	Black	☐	Dove Gray
☐	Cappuccino	☐	Greystone