

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

AM-BAX

PELVIC HEALTH TREATMENT TABLE

STANDARD FEATURES:

- Load capacity: 500 lbs
- 27"W x 72"L x 19-38"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1" High density foam
- Rounded corners
- (4) Total locking casters
- Steel base with dual foot bar control
- Electric elevating midsection
- Adjustable armrests
- Head section angle ranges from +25° to -45°
- Lower extremity angle ranges from 0 to 75°
- White powder coat base



AM-BAX5600


TABLE CONFIGURATIONS:

☐ AM-BAX5600 - 5 Section with Electric Elevating Midsection and Adjustable Armrests 27"W x 72"L x 19-38"H
QTY:



OPTIONS:

☐ AM-27 - Face Plug
☐ AM-91 - Add 1" Extra Soft Foam
☐ 19855 - Add-on Belt & Brackets
☐ SPCO - Special Vinyl Color (Provided by Customer)
☐ SPBASE - Special Base Color (Select Custom Base Color)



STANDARD BASE COLOR:

☐	Office White
-----------------------	--------------

CUSTOM BASE COLORS:

☐	Black
☐	Steel Gray
☐	Metallic Silver

STANDARD UPHOLSTERY CHOICES:

☐	Imperial Blue	☐	Patina
☐	Blue Ridge	☐	Taupe
☐	Forest Green	☐	Burgundy
☐	Black	☐	Dove Gray
☐	Cappuccino	☐	Greystone



800.701.5122
armedicamfg.com