

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

AM-BA SERIES

34" BARIATRIC HI-LO TREATMENT TABLES

STANDARD FEATURES:

- Load capacity: 500 lbs
- 34"W x 76"L x 18-37"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1.5" High density foam
- Rounded corners
- Steel base with dual foot bar control
- White powder coat base
- (4) Total locking casters
- ADA Compliant



AM-BA234

TABLE CONFIGURATIONS:

<input type="radio"/> AM-BA134 - 1 Section	<input type="radio"/> AM-BA234 - 2 Section	<input type="radio"/> AM-BA334 - 3 Section
QTY	QTY	QTY



TABLE FEATURES:

<p>Flat Cushion</p>	<p>Lower Extremity Angle Ranges from 0 to 78°</p>	<p>Lower Extremity Angle Ranges from 0 to 78°</p> <p>Head Section Angle Ranges from +45° to -90°</p>
AM-BA134	AM-BA234	AM-BA334



800.701.5122
armedicamfg.com

OPTIONS:

<input type="radio"/> AM-17 - Face Plug
<input type="radio"/> AM-91 - Add 1" Extra Soft Foam
<input type="radio"/> 19859 - Add-on Belt & Brackets
<input type="radio"/> SPCO - Special Vinyl Color (Provided by Customer)
<input type="radio"/> SPBASE - Special Base Color (Select Custom Base Color)
<input type="radio"/> 12088 - Side Panel 9 x 33.25 Perforated
<input type="radio"/> 12075 - Side Panel 9 x 33.25 Flat

STANDARD UPHOLSTERY CHOICES:

<input type="radio"/> Imperial Blue	<input type="radio"/> Taupe
<input type="radio"/> Blue Ridge	<input type="radio"/> Burgundy
<input type="radio"/> Forest Green	<input type="radio"/> Dove Gray
<input type="radio"/> Black	<input type="radio"/> Greystone
<input type="radio"/> Cappuccino	<input type="radio"/> Merlot
<input type="radio"/> Patina	

STANDARD BASE COLOR:

<input type="radio"/>	Office White
-----------------------	--------------

CUSTOM BASE COLORS:

<input type="radio"/>	Black
<input type="radio"/>	Steel Gray
<input type="radio"/>	Metallic Silver