

Name:	
School/Organization:	
Dealer/Rep:	

Phone:	
E-mail:	
Sales Rep:	Quote:

# CHANGING TABLES

## STANDARD FEATURES:

- Load capacity: 500 lbs
- 34 oz. Heavy-duty vinyl with Permablok 3<sup>®</sup> bacterial protection
- 1.5" High density foam
- Rounded corners
- Single foot control pedal
- White powder coat base
- (4) Total locking casters

## TABLE CONFIGURATIONS:

<input type="radio"/> <b>AM-SX1060</b> - AM-SX 25"W 1 Section Top with Side Rails 25"W x 60"L x 18-37"H	<input type="radio"/> <b>AM-SX1072</b> - AM-SX 25"W 1 Section Top with Side Rails 25"W x 72"L x 18-37"H
QTY	QTY

## TABLE FEATURES:



## STANDARD UPHOLSTERY CHOICES:

<input type="radio"/> Imperial Blue	<input type="radio"/> Taupe
<input type="radio"/> Blue Ridge	<input type="radio"/> Burgundy
<input type="radio"/> Forest Green	<input type="radio"/> Dove Gray
<input type="radio"/> Black	<input type="radio"/> Greystone
<input type="radio"/> Cappuccino	<input type="radio"/> Merlot
<input type="radio"/> Patina	



AM-SX1072

## STANDARD BASE COLOR:

Office White

## CUSTOM BASE COLORS:

<input type="radio"/> Black
<input type="radio"/> Steel Gray
<input type="radio"/> Metallic Silver

## OPTIONS:

<input type="radio"/> <b>AM-804</b> - Paper Holder with Cutter (Accommodates 18-21" Paper Rolls)
<input type="radio"/> <b>AM-SX60</b> - Add Base SX with Casters & 6" Clearance (Standard on AM-SX2572 & AM-SX3076, No Casters) *Adding this option will increase the lowest height to 22"H
<input type="radio"/> <b>AM-91</b> - Add 1" Extra Soft Foam
<input type="radio"/> <b>SPCO</b> - Special Vinyl Color (Provided by Customer)
<input type="radio"/> <b>SPBASE</b> - Special Base Color (Select Custom Base Color)
<input type="radio"/> <b>19855</b> - Add-On Belt and Brackets



800.701.5122  
armedicamfg.com