

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

MAT PLATFORM TABLES

WITH ELECTRIC HEIGHT ADJUSTMENT

Conforms to ANSI/UL Std. IEC 60601,
Certified to CSA Std. C22.1, Std. 60601.1

ETL CLASSIFIED



Intertek
4009555

STANDARD FEATURES:

- Load capacity: 1000 lbs
- Electric elevation: 20-30" H
- ETL Certified
- Indestructible tubing to prevent kinks
- Light up display to indicate standby mode/power applied
- Audible warning with delay safety features when device lowers
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 2" High density foam
- Single mounted hand control
- Rounded corners




AMTM-640

STANDARD CONFIGURATIONS:

☐ AMTM-640 - 48"W x 84"L x 20-30"H	☐ AMTM-641 - Adjustable Backrest (Raises to 80") 48"W x 84"L x 20-30"H	☐ AMTM-650 - 60"W x 84"L x 20-30"H	☐ AMTM-660 - 72"W x 96"L x 20-30"H
QTY:	QTY:	QTY:	QTY:

OPTIONS:

☐ AMTM-03611-FS - Substitute a Foot Switch for Single Mounted Hand Control
QTY:
☐ AMTM-03611-FS-2 - Add a Second Foot Switch (Not available with Hand Control)
QTY:
☐ AMTM-03619-BMS-2 - Add a Second Hand Control (Not available with Foot Switch)
QTY:
☐ AMSBO-047 - Side Bumpers Mat Table, 4' x 7' (Not available with AMTM-641)
QTY:
☐ AMSBO-057 - Side Bumpers Mat Table, 5' x 7'
QTY:
☐ AMSBO-068 - Side Bumpers Mat Table, 6' x 8'
QTY:



AMSBO-047

STANDARD BASE COLOR:

☐	Office White
-----------------------	--------------

CUSTOM BASE COLORS:

☐	Black
☐	Steel Gray
☐	Metallic Silver

STANDARD UPHOLSTERY CHOICES:

☐	Imperial Blue	☐	Taupe
☐	Blue Ridge	☐	Burgundy
☐	Forest Green	☐	Dove Gray
☐	Black	☐	Greystone
☐	Cappuccino	☐	Merlot
☐	Patina		



800.701.5122
armedicamfg.com