

Name:	
School/Organization:	
Dealer/Rep:	

Phone:	
E-mail:	
Sales Rep:	Quote:

PORTABLE PARALLEL BARS

STANDARD FEATURES:


- Load capacity: 600 lbs on the base
 - 300 lbs per bar
- Foldable with casters for easy transport
- Transport mode:
 - 12.25-17.5"W x 60"L x 39.5"H
 - 55° Turn radius
- Extended mode:
 - 45.5"W x 60"L x 28.5-41"H
 - Adjustable handrail: 22.75-26.5"W
- Dual lock mechanism for height adjustment
- (4) 4" Total locking medical grade caster system
- Aluminum handrails
- Welded aluminum frame and walking surface
- Anti-microbial and anti-bacterial powder coat
- Anti-slip platform
- No assembly required




PTPB1000

STANDARD CONFIGURATIONS:

<input type="radio"/>	PTPB1000 - Portable Parallel Bars 12.25 - 17.5"W x 60"L x 39.5"H
QTY:	



Watch Set Up and Take Down
Video, Get User Manual, and
More...



800.701.5122
armedicamfg.com