



## AM-644 Steel Mat Platform

- 1-Section Steel Mat Platform
- 84"L x 48"W
- Fixed Height: 20"
- 2" Firm Density Foam
- 900lbs. Capacity

Also available:  
 AM-654 84"L x 60"W  
 AM-664 96"L x 72"W



## AM-641 Steel Mat Platform

- 2-Section Steel Mat Platform
- 84"L x 48"W
- Min. Height: 20"
- Max Height: 30"
- 2" Firm Density Foam
- Backrest: 29" x 48" with Gas Spring Assist. Locks in Any Position
- Body: 55" x 48"
- 900lbs. Capacity
- Electric Hi-Lo Steel Frame

Also available:  
 AM-643 Adjustable Hi-Lo  
 AM-645 Fixed Height