

AM SERIES

Hi-Lo Treatment Tables



AM-300 Treatment Table

- 3-Section Treatment Table
- 76"L x 27"W
- Min. Height: 18"
- Max Height: 37"
- 1 1/2" Firm Density Foam
- Head: 15" x 27" with contoured face opening
- Center: 22" x 27"
- Foot: 39" x 27"
- 400lbs. Capacity
- Power-assisted Caster System

Clinicians appreciate the full height adjustment and section mobility options so the patient can be placed in the proper position for maximum results with minimum effort.

All tables have these features:

- Comfort and easy access getting on and off of the table
- 34oz. antibacterial vinyl
- 18 month warranty
- Heavy duty welded steel construction