



## AM-400 Traction Table

- |   |   |  |  |
|---|---|--|--|
| <ul style="list-style-type: none"> <li>• 4-Section Traction Table</li> <li>• 76"L x 27"W</li> <li>• 99" Overall Length</li> </ul> | <ul style="list-style-type: none"> <li>• Min. Height: 18"</li> <li>• Max Height: 37"</li> <li>• 1 1/2" Firm Density Foam</li> </ul> | <ul style="list-style-type: none"> <li>• Head: 15" x 27" with contoured face opening</li> <li>• Thoracic: 17.5" x 27"</li> <li>• Center: 10" x 27"</li> <li>• Foot: 33" x 27"</li> </ul> | <ul style="list-style-type: none"> <li>• 400lbs. Capacity</li> <li>• Power-assisted Caster System</li> </ul> |
|---|---|--|--|



## AM-450 Traction Table

- |   |   |  |  |
|---|---|--|--|
| <ul style="list-style-type: none"> <li>• 4-Section Traction Table</li> <li>• 3-Piece Head</li> <li>• 76"L x 27"W</li> <li>• 99" Overall Length</li> </ul> | <ul style="list-style-type: none"> <li>• Min. Height: 18"</li> <li>• Max Height: 37"</li> <li>• 1 1/2" Firm Density Foam</li> </ul> | <ul style="list-style-type: none"> <li>• Head: 15" x 13" with contoured face opening</li> <li>• Arms: 12.5" x 5.5"</li> <li>• Thoracic: 17.5" x 27"</li> <li>• Center: 10" x 27"</li> <li>• Foot: 33" x 27"</li> </ul> | <ul style="list-style-type: none"> <li>• 400lbs. Capacity</li> <li>• Power-assisted Caster System</li> </ul> |
|---|---|--|--|